

Quick Reference Card

Muse® Mitopotential Kit MCH100110

To determine the percentages of cells exhibiting a change in mitochondrial depolarization and cell death

For Research Use Only. Not for use in diagnostic procedures.

Storage Conditions

Store the Muse® MitoPotential Dye at -20°C. Store the Muse MitoPotential 7-AAD Reagent and 1X Assay Buffer at 2 to 8°C.

Kit Components

- Muse® MitoPotential Dye (Part No. 4700-1580, 100 tests/vial)
- Muse MitoPotential 7-AAD (Part No. 4700-1585, 100 tests/vial)
- 1X Assay Buffer (Part No. 4700-1330, 100 tests/vial)

Materials Recommended

- Guava® Muse® Cell Analyzer
- Cell suspension, untreated and treated to undergo apoptosis
- Micropipettors
- Disposable micropipettor tips
- Microcentrifuge tubes with screw caps, 1.5 mL (VWR Catalog No. 16466-030, or equivalent)
- Vortex mixer

Assay Protocol


Culture cells, including positive and negative controls, by desired method.



Prepare Muse® MitoPotential working solution by diluting the dye 1:1000 with 1X Assay Buffer..



Prepare cell samples for incubation with Muse MitoPotential working solution.



The latest version of Muse software, which includes all assay modules, as well as the kit user's guide, can be found at www.luminexcorp.com/flowkits.

Expected Results


Figures A and B show an example of results obtained using the Muse® MitoPotential Kit to stain Jurkat cells treated with staurosporine to induce depolarization.

Events in each of the four quadrants are as follows:


- Lower-left quadrant: live cells with depolarized mitochondrial membrane [MitoPotential (-) and 7-AAD (-)]
- Lower-right quadrant: live cells with intact mitochondrial membrane [MitoPotential (+) and 7-AAD (-)]
- Upper-right quadrant: dead cells with depolarized mitochondrial membrane [MitoPotential (+) and 7-AAD (+)]
- Upper-left quadrant: dead cells with intact mitochondrial membrane [MitoPotential (-) and 7-AAD (+)]

Figures A and B display example data - results obtained with the Muse MitoPotential software module for Jurkat cells stained with the Muse MitoPotential Kit and acquired on the Guava® Muse Cell Analyzer. Figure A shows results without dot plots, while Figure B shows results with optional dot plots. The statistics show the cells/mL in the stained cell sample and the percentages of each population. The first plot in Figure B shows MitoPotential vs. Cell Size and the second plot shows Viability vs. MitoPotential, providing data on four cell populations - Live, Depolarized/Live, Depolarized/Dead, and Dead cells.

Figures A and B



A



B

Related Products

For Research Use Only. Not for use in diagnostic procedures.

Muse® System Check Kit - MCH100101

Muse® Count & Viability Kit - MCH100102

Muse® Count & Viability Kit - MCH600103

Muse® Count & Viability Kit (200X) - MCH100104

Muse® Annexin & Dead Cell Kit - MCH100105

Muse® Cell Cycle Kit - MCH100106

Muse® Cell Dispersal Reagent - MCH100107

Muse® Caspase-3/7 Kit - MCH100108

Muse® MultiCaspase Kit - MCH100109