

Quick Reference Card

Muse® Nitric Oxide Kit MCH100112

For the detection of Nitric Oxide (NO) activity in combination with cell death in cellular samples
For Research Use Only. Not for use in diagnostic procedures.

Storage Conditions

Store the Muse® Nitric Oxide Reagent at -15° to -25°C, protected from light.

Store the Muse 7-AAD and 1X Assay Buffer at 2° to 8°C, protected from light.

Kit Components

- Muse® Nitric Oxide Reagent (Part No. 4700-1666, 100 tests/vial)
- Muse 7-AAD Reagent (Part No. 4700-1673, 100 tests/vial)
- 1X Assay Buffer (Part No. 4700-1330, 100 mL/vial)

Materials Recommended

- Guava® Muse® Cell Analyzer
- Cell suspension; untreated and treated
- Micropipettors
- Disposable micropipettor tips
- Microcentrifuge tubes with screw caps, 1.5 mL (VWR Catalog No. 16466-030, or equivalent)
- Muse Count & Viability CDR (Catalog No. MCH100107), optional
- Vortex mixer

Assay Protocol

Culture cells, including positive and negative controls by desired method.



Prepare cell samples in 1X Assay Buffer at 1×10^6 to 1×10^7 cells/mL for incubation with Muse® Nitric Oxide working solution.



Dilute Muse Nitric Oxide Reagent 1:1000 with 1X Assay Buffer to make Muse® Nitric Oxide working solution.



Add 2 µL of Muse® 7-AAD to 88 µL of 1X Assay Buffer to make Muse 7-AAD working solution.



Add 100 µL of Muse® Nitric Oxide working solution to 10 µL of cells.

Add 90 µL of Muse® 7-AAD working solution.



Incubate at 37°C for 30 minutes.



Mix thoroughly.

Run on Muse® Cell Analyzer.



NOTE: A detailed kit user's guide can be found at www.luminexcorp.com/flowkits (search by Catalog No. MCH100112).

Expected Results


The figures below show an example of results using the Muse® Nitric Oxide Kit to stain RAW 264.7 cells treated with Lipopolysaccharide and mouse IFN γ to induce nitric oxide (NO) expression.

Events in each of the four quadrants are as follows:

- LL quadrant: Viable Cells with no Nitric Oxide Activity, [Nitric Oxide(-) and Dead Cell Marker (-)]
- LR quadrant: Live Cells with Nitric Oxide Activity, [Nitric Oxide(+) and Dead Cell Marker (-)]
- UR quadrant: Dead Cells with Nitric Oxide Activity, [Nitric Oxide(+) and Dead Cell Marker (+)]
- UL quadrant: Dead Cells with no Nitric Oxide Activity [Nitric Oxide(-) and Dead Cell Marker (+)]

Figures A and B. Results obtained from RAW264.7 cells treated with 100 nM LPS and 100 U/mL mouse IFN γ for 22 hours and stained with the Muse Nitric Oxide Kit, then acquired on the Guava® Muse Cell Analyzer.

Figure A shows results without dot plots, while Figure B shows the same results with optional dot plots. The statistics show the cells/mL in the stained cell sample and the percentages of each population. The first plot in B shows Nitric Oxide vs. Cell Size and the second plot shows Nitric Oxide vs. Viability.



For more information, refer to the comprehensive user's guide, which can be found at www.luminexcorp.com/flowkits (search for Catalog Code, MCH100112).

Related Products

For Research Use Only. Not for use in diagnostic procedure.

- Muse® System Check Kit - MCH100101
- Muse® Count & Viability Kit (40 mL) - MCH100102
- Muse® Annexin V & Dead Cell Kit - MCH100105
- Muse® Caspase-3/7 Kit - MCH100108
- Muse® MultiCaspase Kit - MCH100109
- Muse® MitoPotential Kit - MCH100110
- Muse® Oxidative Stress Kit - MCH100111
- Muse® Ki67 Proliferation Kit - MCH100114

The latest version of Muse software, which includes all assay modules, as well as the kit user's guide, can be found at www.luminexcorp.com/flowkits.

Guava and Muse are trademarks of Luminex Corporation, registered in the U.S. and other countries.

Part No. 4600-3453 Rev C, February 2020, Printed in the USA.

© 2020 Luminex Corporation. All rights reserved.