
BeadBug™-6

Microtube Homogenizer

Benchmark
Scientific

PO Box 709 Edison, NJ 08817 USA

Ph: 908 769-5555 EM: info@benchmarkscientific.com

Foreword

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Prefilled 2.0ml tubes

D1031-01	Prefilled 2.0ml tubes, Silica (Glass) Beads, 0.1mm Acid Washed, 50pk
D1031-05	Prefilled 2.0ml tubes, Silica (Glass) Beads, 0.5mm Acid Washed, 50pk
D1031-10	Prefilled 2.0ml tubes, Silica (Glass) Beads, 1.0mm Acid Washed, 50pk
D1032-01	Prefilled 2.0ml tubes, Zirconium Beads, 0.1mm Triple-Pure - High Impact, 50pk
D1032-05	Prefilled 2.0ml tubes, Zirconium Beads, 0.5mm Triple-Pure - High Impact, 50pk
D1032-10	Prefilled 2.0ml tubes, Zirconium Beads, 1.0mm Triple-Pure - High Impact, 50pk
D1032-15	Prefilled 2.0ml tubes, Zirconium Beads, 1.5mm Triple-Pure - High Impact, 50pk
D1032-30	Prefilled 2.0ml tubes, Zirconium Beads, 3.0mm Triple-Pure - High Impact, 50pk
D1032-SK	Prefilled - Triple-Pure Starter Kit, 10 each of 0.1, 0.5, 1.0, 1.5 and 3.0mm
D1033-28	Prefilled 2.0ml tubes, Stainless Steel, 2.8mm Acid Washed, 50pk
D1033-30G	Prefilled 2.0ml tubes, Garnet Shards and one 6 mm Zirconium Bead, 50pk
D1034-MX	Prefilled 2.0ml tubes, 0.1mm Silica, 1.4 mm Zirconium & 4mm Silica Beads, 50pk
D1032-60	Prefilled 2.0ml tubes, 6 mm Ceria Based Zirconium Oxide Satellite, 50pk
D1034-28	Prefilled 5.0ml tubes, Stainless Steel beads, 2.8mm Acid Washed, 50pk**

Empty Tubes

D1031-T20	2.0ml tubes (empty) pack of 1000 with caps and sealing ring
D1031-T21	2.0ml tubes (empty) pack of 50 with caps and sealing ring
D1031-RF	2.0ml reinforced tubes (empty) pack of 500 with cap and sealing ring
C1005-SC5	5.0ml tubes (empty), pack of 500 with cap**

Bulk Beads

D1131-01	Bulk Beads, Silica (glass), 0.1mm, acid washed, 200g
D1131-05	Bulk Beads, Silica (glass), 0.5mm, acid washed, 200g
D1131-10	Bulk Beads, Silica (glass), 1.0mm acid washed, 200g
D1132-01TP	Bulk Beads, Zirconium, 0.1mm, Triple-Pure Molecular Biology Grade, 250g
D1132-05TP	Bulk Beads, Zirconium, 0.5mm, Triple-Pure Molecular Biology Grade, 250g
D1132-10TP	Bulk Beads, Zirconium, 1.0mm, Triple-Pure Molecular Biology Grade, 250g
D1132-15TP	Bulk Beads, Zirconium, 1.5mm, Triple-Pure Molecular Biology Grade, 250g
D1132-30TP	Bulk Beads, Zirconium, 3.0mm, Triple-Pure Molecular Biology Grade, 300g
D1133-28	Bulk Beads, Stainless Steel, 2.8mm, acid washed, 1,000/pk
D1132-60	Bulk Beads, 6 mm Zirconium Oxide, Ceria Stabilized, 50/pk
D1133-G	Garnet shards, bulk, acid washed, 250g bottle

5ml Tube Holder

D1036-A5	Holder for 1x5ml tubes in BeadBug-6, pack of 2
----------	--

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Custom fills available for quantity orders.

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	Bead material	Silica			Zirconium					Steel	Garnet with Zr Satellite	Silica and Zirconium
	Bead diameter	0.1mm	0.5mm	1.0mm	0.1mm	0.5mm	1.0mm	1.0mm	3.0mm	2.8mm	0.5/6mm	Mixed
Animal	Adipose							X				
	Artery							X	X			
	Bone									X		
	Brain							X	X			
	Insect - fly								X	X		
	Insect - Tick								X	X		
	Insect - Mosquito						X	X	X			
	Gonad							X	X			
	Hair								X	X	X	
	Heart							X	X			
	Kidney							X				
	Liver											
	Lung							X	X			
	Muscle								X	X		
	Pancreas							X	X			
	Skin								X		X	
Spleen							X	X				
Teeth								X	X			
Microbial	Algae		X			X						
	Bacteria	X	X		X	X						
	Fungi		X	X		X	X					
	Mold		X			X						
	Spores		X	X		X	X					
	Yeast		X	X		X	X					
Plant and Soil	Leaves								X	X		
	Nuts									X		
	Plant, soft tissue								X	X		
	Plant, tough tissue								X	X	X	
	Seeds								X	X		
	Stems								X	X		
	Soil							X				X
	Coral									X		
Fecal	Stool							X				X
	Scat								X	X		X

Above is a partial list of sample types and suggested beads. For additional information and technical support, contact Benchmark Scientific at info@BenchmarkScientific.com.

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Problems and actions

No.	Problem	Possible cause	Solution(s)
1	No display on the screen	No power	Check the power supply and that the unit is plugged in properly
		On/Off button broken	Change button
		Others	Contact supplier
2	Velocity measurement inaccurate	Sensor broken	Contact supplier
		Electric motor broken	
3	No alarm when the cover open when running	Sensor broken	Contact supplier
4	"E101" displayed with alarming "Du..."	Overspeed alarm	Contact supplier. Press any key to silence the alarm.
5	"E104" displayed with alarming "Du..."	Locked-rotor alarm	Contact supplier. Press any key to silence the alarm
6	"E501" displayed with alarming "Du..."	Cover open when running	Close cover. Press any key to silence the alarm
7	"E503" displayed with alarming "Du..."	Cover didn't close before running	Close the cover. Press any key to silence the alarm.
8	Buttons don't work	Faulty control panel	Contact supplier