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Instruction for Use TPP Tissue Culture Dish

TPP tissue culture dishes are for the cell/tissue cultivation and growth in manual use. The circumferential gripping ring allows a secure grip of the unit (bowl + lid). The growth surface - but not the side wall of the dish - is opto-mechanically activated for optimal cell adhesion and growth.

The tissue culture dish is for single use.

Note:

Follow national regulations when handling biological material, use appropriate protective clothing

Take note of the regulations of aseptic working during the work process.

Instructions:

- Check the packaging material and the product for intactness. Use only faultless material.
- Check the expiry date (EXP) on label and packaging. Only use products with valid EXP.
- Open the packaging in a sterile environment and remove a complete system (lid with base) for use. The circumferential grip ring allows a secure grip of the cell culture dish and thus avoids unintentional lifting of the lid (risk of contamination).
- Remove the lid and fill the dish with the medium and inoculum according to your laboratory routine.
- To achieve an uniform cell growth, the following aspects must be observed:
 - Resuspend the cells well with the medium to avoid cell lumps.

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Advice and Recommendations:

- The stacking rim provides a safe stability of several TPP cell culture dishes stacked on top of each other.
- For a homogenous cultivation of the cells, levelling of the incubator and continuous supply of medium is essential.
- Avoid and/or eliminate vibration sources near the incubator to ensure uniform cell growth.
- Centrifugation of the dish is not recommended.

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Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch.