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Instruction

The PCV tube is intended for the determination of biomass in cell culture suspensions, expressed as packed cell volume.

- Mix the cell suspension thoroughly in TubeSpin® Bioreactor.
- For harvested adherent cells, pipet up and down to thoroughly break cell clumps post trypsinisation.
- Add 200 – 1000 µl of a well-mixed cell suspension to the PCV tube
- Centrifuge for 1 min at 2'500 gx.
 -> Remark: We recommend using a centrifuge with a swinging bucket rotor accommodating 1.5 / 2 ml tubes.
 Alternative: either place the tubes in 15 ml centrifuge tubes with conical bottom and run them in the appropriate rotor or apply a fixed angle rotor instead.
- Read the volume (µl) of the cell pellet directly from the scale of the capillary or with the "easy read" measuring device
- Calculate the PCV:
 $PCV (\%) \times 100$
- Calculate cell count/density (cells/ml)
 perform an initial parallel determination of both cell count (using the hemacytometer) and PCV (using the PCV tube):
 - a. Determine cell number manually with the hemacytometer.
 - b. Draw 1 ml cell suspension from the same sample to determine PCV volume
 - c. Divide the number of cells per ml (from step a) by the volume of the cell pellet in µl (from step b) to get the number of cells per µl cell pellet.
- ➔ Cell pellet cannot be extracted from the capillary.

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Advice and Recommendation

- *The PCV tube provides a fast, easy to use, reproducible and reliable alternative to manual cell counting. The method also is more accurate (margin of error less than 5%). The PCV tubes enable the determination of the packed cell volume (PCV) in a cell suspension, resulting in absolute data correlating with parameters like protein content, cell count, metabolic activity and others.*
- *Samples with high cell densities must not be diluted, simply reduce sample volume. Samples from the same series should be treated accordingly.*
- *For bacteria-measurement the parameters have to be adjusted, for example higher g-force, smaller volume, etc.*
- *PCV tubes are compatible with rotors from most manufacturers.*

**Gebrauchsanleitung TPP - „easy read“ Messinstrument /
 Instruction for Use TPP „easy read“ Measuring Device**

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
Instruction

The “easy read” measuring device is intended to determine the volume of the PCV pellet (PCV = packed cell volume).

1. Put the PCV tube after centrifugation into the “easy read”
2. Move the magnifying lens to the right / left until the upper rim of the pellet is leveled out with the cross lines.
3. Read the volume (µl) of the cell pellet with the “easy read” measuring device

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Technical Data

Materials

Measuring device aluminum
 Magnifying lens glass

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Measurements

Volume graduation in µl
 Length x width x height mm
 Magnifying lens:
 Length x width x height mm

87010

0.0 – 5.0
 320 x 15 x 13.5
 30 x 33 x 38