Challenge

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We do not advocate specific treatments or approaches. We are simply sharing the most recent evidence from the medical community to help other health workers during the SARS-CoV-2 pandemic.

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References

- Miscellaneous inactivating agents - guideline for disinfection and sterilization in healthcare facilities. <https://www.cdc.gov/infectioncontrol/guidelines/disinfection/disinfection-methods/miscellaneous.html> (2016).
- Hijnen, W. A. M., Beerendonk, E. F. & Medema, G. J. Inactivation credit of UV radiation for viruses, bacteria and protozoan (oo)cysts in water: A review. in *Water Research* vol. 40 3–22 (2006).
- Zaffina, S. et al. Accidental Exposure to UV Radiation Produced by Germicidal Lamp: Case Report and Risk Assessment. *Photochem. Photobiol.* 88, 1001–1004 (2012).
- International Commission on Non-Ionizing Radiation Protection. ICNIRP statement-Protection of workers against ultraviolet radiation. *Health Phys.* 99, 66–87 (2010).
- Ridgway, H. F. & Olson, B. H. Chlorine resistance patterns of bacteria from two drinking water distribution systems. *Appl. Environ. Microbiol.* 44, 972–987 (1982).
- Darnell, M. E. R., Subbarao, K., Feinstone, S. M. & Taylor, D. R. Inactivation of the coronavirus that induces severe acute respiratory syndrome, SARS-CoV. *J. Virol. Methods* 121, 85–91 (2004).
- McDevitt, J. J., Milton, D. K., Rudnick, S. N. & First, M. W. Inactivation of Poxviruses by Upper-Room UVC Light in a Simulated Hospital Room Environment. *PLoS ONE* 3, e3186 (2008).
- Park, G. w., Linden, K. g. & Sobsey, M. d. Inactivation of murine norovirus, feline calicivirus and echovirus 12 as surrogates for human norovirus (NoV) and coliphage (F+) MS2 by ultraviolet light (254 nm) and the effect of cell association on UV inactivation. *Lett. Appl. Microbiol.* 52, 162–167 (2011).
- Steinmann, E. et al. Two pathogen reduction technologies—methylene blue plus light and shortwave ultraviolet light effectively inactivate hepatitis C virus in blood products. *Transfusion (Paris)* 53, 1010–1018 (2013).
- Xiong, P. & Hu, J. Inactivation/reactivation of antibiotic-resistant bacteria by a novel UVA/LED/TiO₂ system. *Water Res.* 47, 4547–4555 (2013).
- Zou, S. et al. Inactivation of the novel avian influenza A (H7N9) virus under physical conditions or chemical agents treatment. *Virol. J.* 10, 289 (2013).
- Matsuura, S. & Ishikura, S. Suppression of Tomato mosaic virus disease in tomato plants by deep ultraviolet irradiation using light-emitting diodes. *Lett. Appl. Microbiol.* 59, 457–463 (2014).
- Beck, S. E., Wright, H. B., Hargy, T. M., Larason, T. C. & Linden, K. G. Action spectra for validation of pathogen disinfection in medium-pressure ultraviolet (UV) systems. *Water Res.* 70, 27–37 (2015).
- Lindsley, W. G. et al. Effects of Ultraviolet Germicidal Irradiation (UVGI) on N95 Respirator Filtration Performance and Structural Integrity. *J. Occup. Environ. Hyg.* 12, 509–517 (2015).
- Aubry, M., Richard, V., Green, J., Broult, J. & Musso, D. Inactivation of Zika virus in plasma with amotosalen and ultraviolet A illumination. *Transfusion (Paris)* 56, 33–40 (2016).
- Bae, K. S., Shin, G.-A., Bae, K. S. & Shin, G.-A. Inactivation of various bacteriophages by different ultraviolet technologies: Development of a reliable virus indicator system for water reuse. *Environ. Eng. Res.* 21, 350–354 (2016).
- Bedell, K., Buchaklian, A. & Perlman, S. Efficacy of an automated multi-emitter whole room UV-C disinfection system against Coronaviruses MHV and MERS-CoV. *Infect. Control Hosp. Epidemiol.* 37, 598–599 (2016).
- Cap, A. P. et al. Treatment of blood with a pathogen reduction technology using UV light and riboflavin inactivates Ebola virus in vitro. *Transfusion (Paris)* 56, S6–15 (2016).
- Faddy, H. M. et al. Riboflavin and ultraviolet light: impact on dengue virus infectivity. *Vox Sang.* 111, 235–241 (2016).
- Song, K., Mohseni, M. & Taghipour, F. Application of ultraviolet light-emitting diodes (UV-LEDs) for water disinfection: A review. *Water Res.* 94, 341–349 (2016).

21. Stibich, M. & Stachowiak, J. The microbiological impact of pulsed xenon ultraviolet disinfection on resistant bacteria, bacterial spore and fungi and viruses. *South. Afr. J. Infect. Dis.* 31, 12–15 (2016).
22. Carratalà, A. et al. Experimental adaptation of human echovirus 11 to ultraviolet radiation leads to resistance to disinfection and ribavirin. *Virus Evol.* 3, (2017).
23. Fryk, J. J. et al. Reduction of Zika virus infectivity in platelet concentrates after treatment with ultraviolet C light and in plasma after treatment with methylene blue and visible light. *Transfusion (Paris)* 57, 2677–2682 (2017).
24. Kim, D.-K., Kim, S.-J. & Kang, D.-H. Inactivation modeling of human enteric virus surrogates, MS2, Q β , and Φ X174, in water using UVC-LEDs, a novel disinfecting system. *Food Res. Int.* 91, 115–123 (2017).
25. Maria, F. S. et al. Inactivation of Zika virus in platelet components using amotosalen and ultraviolet A illumination. *Transfusion (Paris)* 57, 2016–2025 (2017).
26. Mills, D., Harnish, D. A., Lawrence, C., Sandoval-Powers, M. & Heimbuch, B. K. Ultraviolet germicidal irradiation of influenza-contaminated N95 filtering facepiece respirators. *Am. J. Infect. Control* 46, e49–e55 (2018).
27. Emig, E. et al. Efficacy of a novel ultraviolet light-emitting diode device for decontamination of shared pens in a health care setting. *Am. J. Infect. Control* 48, 100–102 (2020).
28. Lowe, J. et al. N95 Filtering Facemask Respirator Ultraviolet Germicidal Irradiation (UVGI) Process for Decontamination and Reuse. <https://www.nebraskamed.com/sites/default/files/documents/covid-19/n-95-decon-process.pdf> (2020).
29. Recommended Guidance for Extended Use and Limited Reuse of N95 Filtering Facepiece Respirators in Healthcare Settings. The National Institute for Occupational Safety and Health (NIOSH) - Pandemic Planning <https://www.cdc.gov/niosh/topics/hcwcontrols/recommendedguidanceextuse.html> (2018).
30. Blachere, F. M. et al. Assessment of influenza virus exposure and recovery from contaminated surgical masks and N95 respirators. *J. Virol. Methods* 260, 98–106 (2018).
31. 2019 Novel Coronavirus and COVID-19 Disease Outbreak. Worker Health & Safety - Novel Coronavirus and COVID-19 Outbreak - 3M Personal Protective Equipment (PPE) Considerations Technical Bulletin Revision 18 https://www.3m.com/3M/en_US/worker-health-safety-us/all-stories/full-story-detail/?storyid=8855304f-01cb-4af2-8937-83096cdb4113 (2020).
32. Price, A and Chu LF. COVID-19 - Addressing COVID-19 Face Mask Shortages (v 1.3), Learnly Anesthesia/Stanford AIM Lab COVID-19 Evidence Service, Stanford, California USA (2020).

This document is true and correct at the time of publication; the information within is subject to change. Other documents may supersede this document, including technical modifications and corrections.

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E-mail sales@analytik-jena.co.jp
URL <http://www.analytik-jena.de/jp>

jpn · 04/2020

© Analytik Jena US LLC | Pictures ©: Analytik Jena / freepik
Printout and further use permitted with reference to the source