

2-16P/2-16KL用ローターリスト

アングル (アダプター別売)												
カタログ番号	チューブ数	最大回転数 (RPM)	最大遠心力 (g)	希望販売価格	バケット	入数	希望販売価格	チューブ容量	アダプター	入数	希望販売価格	備考
SGM12104	96	14,000	外側: 21,475 内側: 18,297	¥211,000				0.2mL		不要		
SGM12148	24	15,000 (2-16P) 15,300 (2-16KL)	20,627 (2-16P) 21,460 (2-16PK)	¥161,000				0.2mL	SGM13021	2	¥1,000	
								0.5mL	SGM13002	2	¥1,000	
								1.5mL		不要		
								2.0mL		不要		
SGM12133 (2-16P用)	30	12,000	16,099	¥185,000				0.2mL	SGM13021	2	¥1,000	
								0.5mL	SGM13002	2	¥1,000	
								1.5mL		不要		
								2.0mL		不要		
SGM12132 (2-16KL用)	30	14,000 13,500@4℃	21,913	¥185,000				0.2mL	SGM13021	2	¥1,000	
								0.5mL	SGM13002	2	¥1,000	
								1.5mL		不要		
								2.0mL		不要		
SGM12180	12	14,000 (2-16P) 15,000 (2-16KL)	17,530 (2-16P) 20,124 (2-16KL)	¥161,000				1.5-1.8mL	SGM14080	4	¥7,000	クライオチューブ用 ローターリッド不可
								1.5mL				
								2.0mL	SGM14082	4	¥7,000	
								5.0mL		不要		
SGM12141	10	13,500 (2-16P) 15,300 (2-16KL)	15,485 (2-16P) 19,890 (2-16KL)	¥179,000				10mLラウンド (16.1mmΦ)		不要		
SGM12139	6	13,500 (2-16P) 15,300 (2-16KL)	15,893 (2-16P) 20,414(2-16KL)	¥198,000				10mLラウンド (16.1mmΦ)	SGM13059	2	¥18,000	
								30mLラウンド		不要		
SGM12181	6	10,000 (2-16P) 12,000 (2-16KL)	10,621 (2-16P) 15,294 (2-16KL)	¥227,000				15mLコニカル	SGM13060	2	¥9,000	
								25mL/30mL自立	SGM13064	2	¥17,000	
								50mLラウンド	SGM13079	2	¥11,000	
								50mLコニカル		不要		
SGM12138 (2-16KL用)	4	9,400	9,484	¥229,000				10mLラウンド (16.1mmΦ)	SGM13084	2	¥17,000	2本/アダプター
								17mm	SGM16019	2	¥16,000	
								15mLコニカル	SGM13081	2	¥17,000	
								30mLラウンド (25mmΦ)	SGM13083	2	¥17,000	
								30mLラウンド (23.5mmΦ)	SGM16031	2	¥11,000	
								50mLラウンド	SGM13062	2	¥13,000	
								50mLコニカル	SGM13080	2	¥17,000	
								85mLラウンド		不要		
SGM12071	外側 12 内側 6	4,000	外側: 2,254 内側: 1,842	¥262,000	外側 SGM13011 (標準付属) 内側 SGM13151 (標準付属)	2	¥5,000 ¥55,000	12x75mm	SGM13071	2	¥7,000	
								12x100mm	SGM13072	2	¥7,000	
								16x80mm	SGM13037	2	¥7,000	
								16-17x100mm		不要		
								15mLコニカル		不要		
								25mL/30mL自立	SGM13064	2	¥17,000	
								50mLコニカル		不要		
SGM12072	20	4,000	2,486	¥158,000	SGM13011 (標準付属)	2	¥5,000	12x75mm	SGM13071	2	¥7,000	
								12x100mm	SGM13072	2	¥7,000	
								16x80mm	SGM13073	2	¥7,000	
								16-17x100mm		不要		
								15mLコニカル		不要		
SGM12073	20	4,000	2,486	¥168,000				12x75mm	SGM13071	2	¥7,000	
								12x100mm	SGM13072	2	¥7,000	
								16x80mm	SGM13073	2	¥7,000	
								16-17x100mm		不要		
								15mLコニカル		不要		

