

BMBioリザーバー耐薬品リスト

	20°C	50°C
Acetone	S	S
Ammonium Chloride	E	E
Ammonium Sulfate	E	B
Dimethylsulfoxide	S	S
Ethanol	E	S
Chloroform	S	N
Magnesium Chloride	E	E
Mineral Oil	E	S
Phenol 10%	S	S
Tromethamine	E	E
Edeate Disodium	E	B
Acetic Acid, Glacial	S	S
Benzene	S	S
Boric Acid	E	B
1,3-Butadiene	E	B
Calcium Hydroxide	E	S
Chlorine	S	S
Formaldehyde	E	E
Hexanes	E	E
Hydrochloric Acid	E	S
Hydrogen Peroxide	E	S
Kerosene	E	E
Sulfuric Acid 95%	N	N
Sulfuric Acid 20%	E	E

E = Excellent resistance, no etching

B = Good resistance, little etching after day 28

S = Fair resistance, etching after day 7

N = Not recommended