

		シリカ			ジルコニウム					ステンレス スチール	ガーネット断片 /ジルコニウム	シリカ /ジルコニウム
		0.1 mm	0.5 mm	1.0 mm	0.1 mm	0.5 mm	1.0 mm	1.5 mm	3.0 mm	2.8 mm	0.5/6.0 mm	混合
		D1031-01	D1031-05	D1031-10	D1032-01	D1032-05	D1032-10	D1032-15	D1032-30	D1033-28	D1033-30G	D1034-MX
動物	脂肪							●				
	血管							●	●			
	骨									●		
	脳							●	●			
	性腺							●	●			
	毛髪								●	●	●	
	心臓							●	●			
	腎臓							●	●			
	肝臓							●	●			
	肺							●	●			
	筋肉								●	●		
	膵臓							●	●			
	皮膚								●		●	
	脾臓							●	●			
	歯								●	●		
	ハエ								●	●		
ダニ								●	●			
蚊							●	●	●			
微生物	藻類		●			●						
	細菌	●	●		●	●						
	真菌		●	●		●	●					
	カビ		●			●						
	孢子		●	●		●	●					
	酵母		●	●		●	●					
植物 / 土壌	葉								●	●		
	ナッツ								●	●		
	軟組織								●	●		
	硬組織								●	●	●	
	種子								●	●		
	茎								●	●		
	土壌							●				●
糞便	糞便(ヒト)							●				●
	糞便(動物)								●	●		●