





























































CENTRIFUGE

GYROZEN遠心機とローター

用途に合わせてローターをお選びください。

多目的 高速 遠心機	1248 / 1248R	Rotor Cat.No.	GRF-G-s0.2-64	GRF-G-m2.0-30	GRA-G-15-24A GLB-15A-FA	GRF-G-15-24	GRF-G-50-8	GRF-G-85-6							
		Bucket Cat.No.	-	-	-	-	-	-	-						
		Max. Capacity	64x0.2ml	30x2.0ml	24x15ml	24x15ml	8x50ml	6x85ml							
		Shape													
		Max. RPM/RCF	12,000/9,434(inner), 10,800(outer)	12,000/15,520	4,000/3,134	8,000/8,872	12,000/16,582	10,000/12,298							
		Hole Diameter(mm)	6.5	11.1	18	17.5	29.5	38.5							
	Max. Height for Tube fit(mm)	16	52	125	115	120	122								
	Rotor Cat.No.	GRS-G-100-6						GRS-G-250-4			GRS-G-r250-4				
	Bucket Cat.No.	GLB-b100-100	GLB-b50-100	GLB-10-8-100	GLB-bd15-100	GLB-bdc15-100	GLB-250-250	GLB-c250-250	GLB-60-250	GLB-d20/10-250	GLB-r250-r250	GLP-mw-r250			
	Max. Capacity	6x100ml	6x50ml	6x8x10ml vacutainer	6x2x15ml	6x2x15ml conical	4x250ml	4x250ml conical	4x60ml	4x2x(20&10ml)	4x250ml	4x4 MTPs			
	Shape														
	Max. RPM/RCF	4,000/3,105	4,000/3,032	4,000/3,105	4,000/3,032	4,000/3,105	5,000/5,394	5,000/5,394	5,000/5,394	5,000/4,752(20), 4,472(10)	4,000/3,134	4,000/2,737			
Hole Diameter(mm)	47	30.5	16	17	17	62	60.5	31.5	22.5 (20), 16 (10)	86x70.3	88x128.5				
Max. Height for Tube fit(mm)	130	125	120	120	127	153	153	165	132 (20), 125 (10)	130	80				
多目的 高速 遠心機	1580 / 1580R	Rotor Cat.No.	GRF-G-s0.2-64	GRF-G-m2.0-30	GRA-G-15-24A GLB-15A-FA	GRF-G-15-24	GRF-G-50-8	GRF-G-85-6	GRF-G-250-6						
		Bucket Cat.No.	-	-	-	-	-	-	-						
		Max. Capacity	64x0.2ml	30x2.0ml	24x15ml	24x15ml	8x50ml	6x85ml	6x(250&15ml)						
		Shape													
		Max. RPM/RCF	12,000/9,434(inner), 10,800(outer)	15,000/24,249	4,000/3,134	8,000/8,872	15,000/25,910 12,000/16,582	12,000/17,709 10,000/12,298	8,000 9,588(250), 9,389(15)						
		Hole Diameter(mm)	6.5	11.1	18	17.5	29.5	38.5	62(17)						
	Max. Height for Tube fit(mm)	16	52	130	115	120	122	130(120)							
	Rotor Cat.No.	GRS-G-100-6						GRS-G-250-4			GRS-G-r250-4				
	Bucket Cat.No.	GLB-b100-100	GLB-b50-100	GLB-10-8-100	GLB-bd15-100	GLB-bdc15-100	GLB-250-250	GLB-c250-250	GLB-60-250	GLB-d20/10-250	GLB-r250-r250	GLP-mw-r250			
	Max. Capacity	6x100ml	6x50ml	6x8x10ml vacutainer	6x2x15ml	6x2x15ml conical	4x250ml	4x250ml conical	4x60ml	4x2x(20&10ml)	4x250ml	4x4 MTPs			
	Shape														
	Max. RPM/RCF	4,000/3,105	4,000/3,032	4,000/3,105	4,000/3,032	4,000/3,105	5,000/5,394	5,000/5,394	5,000/5,394	5,000/4,752(20), 4,472(10)	4,000/3,134	4,000/2,737			
Hole Diameter(mm)	47	30.5	16	17	17	62	60.5	31.5	22.5 (20), 16 (10)	86x70.3	88x128.5				
Max. Height for Tube fit(mm)	130	125	120	120	127	153	153	165	132 (20), 125 (10)	130	80				
Rotor Cat.No.							GRS-G-r750-4			GRS-G-750-4					
Bucket Cat.No.							GLB-r750-r750	GLB-500-r750	GLP-mw-r750	GLB-750-750	GLB-b750-750	GLB-500-750	GLP-mw-750		
Max. Capacity							4x500ml conical	400x500ml	4x4 MTPs	4x750ml	4x750ml	4x500ml	4x4 MTPs		
Shape															
Max. RPM/RCF							3,600/2,898	3,600/3,159	3,600/2,485	4,000/3,667		4,000/3,900	4,000/3,068		
Hole Diameter(mm)							99x99	69.2	88x128.5	99	99	69.2	88x128.5		
Max. Height for Tube fit(mm)							135	165	80	142	130	165	80		

小型 遠心機	MINI	Rotor Cat.No.	GRA-s0.2-32	GRF-m1.8-10	GRF-m2.0-12							
		Max. Capacity	32x0.2ml	10x1.8ml cryotube	12x2.0ml							
		Shape										
		Max. RPM/RCF	6,000/1,554(1row), 1,650(2row)	13,500/12,300	13,500/12,300							
		Hole Diameter(mm)	6.5	12.8	11.1							
		Max. Height for Tube fit(mm)	-	50	42							
	1524	Rotor Cat.No.	GRF-s0.2-32	GRF-m2.0-18	GRF-m2.0-24							
		Max. Capacity	32x0.2ml	18x2.0ml filter tube	24x2.0ml							
		Shape										
		Max. RPM/RCF	15,000/14,187	15,000/18,891	15,000/21,206							
		Hole Diameter(mm)	6.5	11.1	11.1							
		Max. Height for Tube fit(mm)	16	65	52							
1536	Rotor Cat.No.	GRF-s0.2-32	GRF-m2.0-24A	GRF-m2.0-36								
	Max. Capacity	32x0.2ml	24x2.0ml	36x2.0ml								
	Shape											
	Max. RPM/RCF	15,000/14,187	15,000/21,206	15,000/21,206								
	Hole Diameter(mm)	6.5	11.1	11.1								
	Max. Height for Tube fit(mm)	16	56	52								
1730R	Rotor Cat.No.	GRF-M-s0.2-64	GRF-M-m2.0-18	GRF-M-m2.0-24	GRF-M-m2.0-30	GRD-M-m2.0-8	GRD-M-m2.0-16	GRF-M-m5.0-12				
	Max. Capacity	64x0.2ml	18x2.0ml filter tube	24x2.0ml	30x2.0ml	8x2.0ml	16x2.0ml	12x5.0ml				
	Shape											
	Max. RPM/RCF	12,500/10,237(1row), 11,722(outer)	15,000/18,891	17,000/27,237	14,000/21,124	13,000/16,343	13,000/15,682	15,000/21,583				
	Hole Diameter(mm)	6.5	11.1	11.1	11.1	11.1	11.1	17				
	Max. Height for Tube fit(mm)	16	65	52	52	50	43	70				

予告なく仕様変更になる場合がございます。

SAMPLE PREPARATION & PROCESSING

SHAKING & INCUBATION

LIQUID HANDLING

CENTRIFUGATION

GEL DOCUMENTATION

GENE AMPLIFICATION

GENE EXPRESSION

CELL BIOLOGY

IMAGING

