

## SoftFit~L ピペットチップ適合表

### ノンフィルターチップ

Pipet-Lite XLS + シングルチャンネルピペット\*  
E4 XLS + 電動シングルピペット\*

Pipet-Lite XLS + マルチチャンネルピペット\*  
E4 XLS + 電動マルチピペット\*  
E4 XLS + アジャスタブルスパーサー\*

チップタイプ	型式	0.1-2 µL	0.5-10 µL	2-20 µL	10-100 µL	20-200 µL	20-300 µL	100-1000 µL	1-10 µL	2-20 µL	5-50 µL	20-200 µL	20-300 µL	100-1200 µL
SoftFit~L 20µL	3722-**	●	●	●					●	●				
	3721-**	●	●	●					●	●				
	3722-05-**	●	●	●					●	●				
	3721-05-**	●	●	●					●	●				
SoftFit~L 200µL	3752-**				●	●					●	●		
	3751-**				●	●					●	●		
	3752-05-**				●	●					●	●		
	3751-05-**				●	●					●	●		
SoftFit~L 300µL	3732-**							●						●
	3731-**							●						●
	3732-05-**							●						●
	3731-05-**							●						●
SoftFit~L 1000µL	3782-**							●						
	3781-**							●						
	3782-05-**							●						
	3781-05-**							●						
SoftFit~L 1200µL	3742-**													●
	3741-**													●
	3742-05-**													●
	3741-05-**													●

### フィルターチップ

チップタイプ	型式	0.1-2 µL	0.5-10 µL	2-20 µL	10-100 µL	20-200 µL	20-300 µL	100-1000 µL	1-10 µL	2-20 µL	5-50 µL	20-200 µL	20-300 µL	100-1200 µL
SoftFit~L 20µL	2749-**	●	●	●					●	●				
	2749-05-**	●	●	●					●	●				
SoftFit~L 200µL	2769-**				●	●					●	●		
	2769-05-**				●	●					●	●		
SoftFit~L 300µL	2739-**						●							●
	2739-0-**						●							●
SoftFit~L 1000µL	2779-**							●						
	2779-05-**							●						
SoftFit~L 1200µL	2789-**													●
	2789-05-**													●

\* ピペットはメトラー・トロード社製品です。

\*\* ラック(HR)もしくは詰め替え(RI)が入ります。